

		Nu-Instruments Plasma MC-ICPMS
EPOND Part-No.	OEM Part-No.	Description
<b>Nebulizers/Spray chambers/accessories</b>		
<b>Quartz Nebulizers</b>		
EP 2050-SILQ-50-05		MEINHARD® SiQ+ Ultra high purity quartz nebulizer
EP 2050-TPQ50-A05		MEINHARD® TQ+ quartz nebulizer
EP 2050-MQ30-A05		MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake
EP 2050-MQ30-A1		MEINHARD® Type A concentric quartz nebulizer, 1 ml/min uptake
<b>Glass Nebulizers</b>		
EP 2050-TRP050-05		MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A05		MEINHARD® Type A glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A1		MEINHARD® Type A glass nebulizer, 1 ml/min uptake
<b>PFA Nebulizer</b>		
ES-2040-28		MicroFlow PFA-ST self-aspirating nebulizer
<b>Poly-Pro Nebulizer</b>		
EP 2550-4040-28		Poly-Pro-ST self-aspirating nebulizer
<b>Nebulizer Cleaner</b>		
PRL02		Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
<b>PFA Cyclonic Spray Chamber</b>		
EP 3193-3111-21		PFA cyclonic spray chamber for Nu AttoM, baffled
<b>Torches/shields</b>		
<b>Quartz Torch Bonnet</b>		
EP 2045-0105	319-311	Quartz bonnet for one-piece torches
<b>Quartz One-Piece Torches with Ball Joint</b>		
EP 2045-0050	319-308	Quartz torch, axial, fixed 1.5mm injector
EP 2045-0050R		Quartz torch, axial, fixed 1.5mm injector(rebuild of ML170005
EP 2045-A050	3201192	Quartz torch, axial, fixed 1.5mm injector
<b>Nu AttoM/Plasma II Demountable Quartz Torch</b>		
EP 1602-1000-93		O-Ring-Free Demountable Precision Quartz Torch for Nu AttoM/Nu Plasma II
<b>Injectors and accessories</b>		
<b>Nu AttoM/Plasma II Demountable injectors</b>		
EP 1233-0150-93		Quartz injector for Nu AttoM, 1.5mm I.D., Nu Plasma II, o-ring-free
EP 1233-0200-93		Quartz injector for Nu AttoM, 2mm I.D., Nu Plasma II, o-ring-free
EP 1203-0100-93		Sapphire injector for Nu AttoM, 1mm I.D., Nu Plasma II, o-ring-free
EP 1203-0150-93		Sapphire injector for Nu AttoM, 1.5mm I.D., Nu Plasma II, o-ring-free
EP 1203-0180-93		Sapphire injector for Nu Attom, 1.8mm I.D., Nu Plasma II, o-ring-free
EP 1213-0100-93		Platinum injector for Nu Attom, 1mm I.D., Nu Plasma II, o-ring-free
EP 1213-0150-93		Platinum injector for Nu Attom, 1.5mm I.D., Nu Plasma II, o-ring-free
EP 1213-0200-93		Platinum injector for Nu Attom, 2mm I.D., Nu Plasma II, o-ring-free
<b>O-Ring-Free Sapphire Injectors, Demountable</b>		
EP 1203-0100-82		Sapphire injector o-ring-free, HF resistant, 1mm
EP 1203-0180-82		Sapphire injector o-ring-free, HF resistant, 1.8mm
<b>O-Ring-Free Platinum Injectors, Demountable</b>		
EP 1213-0180-82		Platinum injector o-ring-free, HF resistant, 1.8mm
EP 1213-0200-82		Platinum injector o-ring-free, HF resistant, 2mm

**Cones****Sampler Cones**

EP 2045-0850	319-285	Ni sampler Cone, Ni base(Version A)
EP 2045-0410	319-541	Ni sampler Cone, Ni bas(Version B)
Nu 1001-GG	3004382	Gasket, Graphite for sampler(Qty 3)
Nu1001A-Ni	319-285	Nickel sampler cone, original style
Nu100A4-Ni	319-497	Nickel sampler cone, WA6, original style
Nu1001B-Ni	319-541	Nickel sampler cone, Version B, 0.9mm orifice
Nu1004B-Ni	319-540	Nickel sampler cone, Version B

**Skimmer Cones**

EP 2045-0840	319-284	Ni Skimmer Cone, Ni base(Original)
EP 2045-0970	319-497	Ni Skimmer Cone, Ni base(Version A)
EP 2045-0540	319-540	Ni Skimmer Cone, Ni base(version B)

**Cone Cleaning Kit & Detergent**

EPMNTKIT1		Cone Cleaning Kit
EPCitranox		Citranox Detergent

**Accessories****Magnetic Sector Multipliers**

EPMAF143		ETP Magnetic Sector Multiplier for Nu Plasma
EPMAF144		ETP Magnetic Sector Multiplier for Nu Plasma AttoM, with filter

**RS12/5 Socket to Torch Base PFA Adapter**

EP 5509-0000		PFA adapter, RS12/5 socket to ESI o-ring-free torch base
EP 5510-0000		PFA adapter, ESI spray chamber to RS12/5 torch base