

**Nebulizers/Spray chambers/accessories**
**Quartz Nebulizers**

EP 2050-SILQ-50-05	MEINHARD® SILQ+Ultra high purity quartz nebulizer
EP 2050-TQP-50-A05	MEINHARD® TQP+ Ultra high purity quartz nebulizer, 0.5 ml/min uptake, 50 PSI
EP 2050-MQ30-A05	MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 30 PSI
EP 2050-MQ30-A1	MEINHARD® Type A quartz nebulizer, 1 ml/min uptake 30 PSI
EP 2050-MQ30-A3	MEINHARD® Type A concentric quartz nebulizer, 3 ml/min uptake
EP 2050-MQ50-A05	MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 50 PSI

**Glass Nebulizers**

EP 2050-TRP50-05	MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A05	MEINHARD® Type A glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A1	MEINHARD® Type A glass nebulizer, 1 ml/min uptake, 30 PSI
EP 2050-M030-A3	MEINHARD® Type A glass nebulizer, 3ml/min uptake
EP 2050-M030-C3	MEINHARD® Type C glass nebulizer, 3ml/min uptake

**PFA Nebulizers**

ES-2040-28	MicroFlow PFA-ST self-aspirating nebulizer
------------	--

**Pylo-Pro Nebulizer**

EP 2550-4040-28	Poly-Pro self-aspirating nebulizer
-----------------	------------------------------------

**PEEK Nebulizer**

EP 3150-0720	PEEK MiraMist nebulizer
--------------	-------------------------

**PFA Sample Capillaries for ST Nebulizers**

EP 2550-2045	Self-aspirating sample capillary for ST Nebulizer, 0.15mm ID(red marker), 20µL/min
EP 2550-2043	Self-aspirating sample capillary for ST Nebulizer, 0.20mm ID(purple marker), 50µL/min
EP 2550-2042	Self-aspirating sample capillary for ST Nebulizers, 0.25mm ID(green marker), 100µL/min
EP 2550-2047	Self-aspirating sample capillary for ST Nebulizers, 0.30mm ID(yellow marker), 200µL/min
EP 2550-2041	Self-aspirating sample capillary for St Nebulizers, 0.50mm ID(orange marker), 400µL/min
EP 2550-2044	Self-aspirating sample capillary for ST Nebulizers, 0.80mm ID(blue marker), 700µL/min
EP 2550-2049	Self-aspirating sample capillary for ST Nebulizers, 1.00mm ID(gray marker), 1mL/min

**Nebulizer Cleaner**

PRL02	Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
-------	--

**pergo Argon Nebulizer Gas Humidifier**

PRG-01	<i>pergo</i> Argon high solids without nebulizer, includes gas line connector kit
PRG-01-MaPlus	<i>pergo</i> Argon high solids including MEINHARD® Type A Plus high performance nebulizer

**Scott Spray Chamber**

EP 2045-0160	3200841	Quartz Scott spray chamber, jacketed for coolant
--------------	---------	--

**Conical Spray Chambers**

EP 2045-0220	3204313	Pyrex conical spray chamber, impact, jacketed
EP 2045-0230	3600374	Quartz conical spray chamber, impact , jacketed
EP 2045-0810	3600170	Quartz conical spray chamber, Peltier cooled
EP 2045-150	4600508	Quartz conical spray chamber, peltier cooled, small drain

**Cyclonic Spray Chamber**

EP 2045-0270	3600093	Quartz cyclonic spray chamber for VG
--------------	---------	--------------------------------------

**Torches/shields****Semi-Demountable Torch**

EP6608-1102 Semi-demountable o-ring-free torch with sapphire injector, HF resistant

**Quartz Torches, Demountable Injector**

EP 2045-0510 3601146 O-ring-free torch with pin, demountable

EP 2045-0100 3600970 Quartz torch, demountable, axial

EP 1202-0301 ESI o-ring-free precision quartz torch, demountable, axial

**Torch Bonnet**

EP 2045-0250 3601149 Torch bonnet for XSERIES 2

**Injectors and accessories****O-Ring-Free Platinum Injectors**

EP 1213-3100 Platinum injector, 1 mm, o-ring-free

EP 1213-3150 Platinum injector, 1.5mm, o-ring-free

EP 1213-3200 Platinum injector, 2mm, o-ring-free

**O-Ring-Free Sapphire Injectors**

EP 1203-3150 Sapphire injector, 1.5mm, o-ring-free, HF-resistant

EP 1203-3180 Sapphire injector, 1.8mm, o-ring-free, HF-resistant

**O-Ring-Free Quartz Injectors**

EP 1233-3100 Quartz injector, 1mm, o-ring-free

EP 1233-3150 Quartz injector, 1.5mm, o-ring-free

EP 1233-3200 Quartz injector, 2mm, o-ring-free

EP 1233-3250 Quartz injector, 2.5mm, o-ring-free

**Cones****Cone Cleaning Kit&Detergent**

EPMNTKIT1 Cone Cleaning Kit

EPCitranox Citranox Detergent

**Nickel Cones**

VG1022 3600811 Ni skimmer cone for Thermo X SERIES 2

VG1021 3600812 Ni sampler cone for Thermo X SERIES 2

**Accessories****Carbon Fiber Supported Autosampler Probes**

ES-5037-3250-080 0.25mm I.D. Carbon fiber support(125mm), 100µl/min uptake, 80cm capillary

ES-5037-3500-080 0.5mm I.D. Carbon fiber support(125mm), 400µl/min uptake, 80cm capillary

ES-5037-3250-100 0.25mm I.D. Carbon fiber support(150mm), 100µl/min uptake, 100cm capillary

ES-5037-3500-100 0.5mm I.D. Carbon fiber support(150mm), 400µl/min uptake, 100cm capillary

**HF Resistant Sample Intro Kits**

EP 2085-SIK1 Complete HF-resistant sample introduction kit

EP 2085-SIK2 Complete high purity sample introduction kit

**Quartz Couplings**

EP 2045-0320 3200727 Quartz coupling, ball/socket, 90° for Scott spray chamber

EP 2045-0300 3600169 Quartz coupling, ball/socket, 17°

**Internal Liquid Fitting**

EP 2080-0500QD Internal liquid fitting for spray chamber drain

**Plasma Screen**

EP 2045-9900 3601219 Torch bonnet for XSERIES 2 torch

**Pipette Tips**

EP 7000-0101 N0777394 100µl PFA PureTip pipette tips, pre-cleaned, pack of 10

EP 7000-1001 N0777395 1000µl PFA PureTip pipette tips, pre-cleaned, pack of 10

**PC<sup>3</sup> Peltier Cooler**

EP 4190-60-25 PC<sup>3</sup> Peltier Cooler, quartz cyclonic spray chamber+PFA-ST nebulizer included

EP 9190-60-29		<b>PC<sup>3</sup>-LT Peltier Cooler</b> PC <sup>3</sup> -LT Peltier Cooler -20°C, quartz cyclonic spray chamber+PFA-ST nebulizer included
EP 2045-0473		<b>Copper Load Coil</b> Copper load coil for Thermo XSERIES 2
VG1023-Cu	3600813	<b>Copper Skimmer Adapter</b> Copper skimmer adapter with screws for Thermo XSERIES 2
EPM44721	1244721	<b>Electron Multiplier</b> MasCom Electron Multiplier for Thermo XSERIES 2